



SOLUTIONS TO HELP



5 REASONS TO CHOOSE EEGDIGITRACK!

1. The most **user-friendly EEG solution**
2. **Innovative** and constantly developed system
3. Reliable and **intuitive**
4. **Top quality** products
5. Efficient and fast **service**, also remote

ELMIKO
LONG-TERM aEEG/CFM
NIRS/SPO₂ MONITORING



FIND US: ELMIKO.PL

aEEG/CFM ELMIKO



CEREBRAL FUNCTION MONITORING
WITH SpO₂ AND rSO₂ MEASUREMENTS



SOLUTIONS TO HELP

FIND US: ELMIKO.PL

aEEG/CFM HEADBOX THE MOST IMPORTANT FEATURES

- / unique function of **automatic color marking** of the record depending on the recorded voltage
- / ability to define your own patterns of activity
- / **advanced** trend and EEG signal measurement ruler
- / synchronized analysis video monitoring from up to two cameras
- / many analyzes available, including: **FFT, CSA/DSA, BSI, 3D mapping, TPM**
- / modern EEG CFM headbox with **EEG, ECG, SpO₂** and rSO₂ measurements
- / visualisation in standard mode or with shading to emphasize values of the EEG trend
- / possibility to define many different alarms
- / integrated amplitude and frequency analysis
- / possibility to perform **bipolar analysis** with use of montages

aEEG/CFM EEGDIGITRACK SOFTWARE HIGHLIGHTS

- / **multi-station** work on a distributed database
- / **user friendly** and convenient interface
- / advanced measuring **ruler**
- / adjustable workplace
- / analysis: **FFT, mapping, TPM, DSA / CSA, BSI** and others
- / meets the requirements of the **GDPR**
- / continuous development of the software
- / **backward compatibility**



ELMIKO BIOSIGNALS sp. z o. o.

📍 Jeżewskiego 5c/7
02-796 Warsaw, Poland

☎ +48 22 644 37 37
☎ +48 600 37 37 90

✉ elmiko@elmiko.pl
elmiko.pl



ISO13485 QUALITY MANAGEMENT SYSTEM AND CE MARK CERTIFIED BY TÜV RHEINLAND®

CREATE TO HELP