


SOLUTIONS TO HELP



5 REASONS TO CHOOSE EEGDIGITRACK!

1. The most **user-friendly EEG** solution
2. **Innovative** and constantly developed system
3. Reliable and **intuitive**
4. **Top quality** products
5. Efficient and fast **service**, also remote

ELMIKO
ROUTINE EEG
ELMIKO EXG-32



FIND US: ELMIKO.PL

EEG DigiTrack

ELMIKO EXG-32



32-CHANNEL EEG

THE PERFECT SOLUTION
FOR ROUTINE EEG TESTS

SOLUTIONS TO HELP

FIND US: ELMIKO.PL





ELMIKO EXG-32 HEADBOX THE MOST IMPORTANT FEATURES

- / 22 EEG reference channels
- / 10 bipolar channels
- / damage-resistant construction
- / **wireless** and USB connectivity
- / possibility to connect **photic stimulator**, **photo detector** and **patient button**
- / **SpO₂** pulse oximeter input
- / signaling of channel impedance

ELMIKO EEGDIGITRACK SOFTWARE HIGHLIGHTS

- / **multi-station** work on a distributed database
- / **user friendly** and convenient interface
- / advanced measuring **ruler**
- / adjustable workplace
- / possible expansion with modules:
Video, PSG, Holter, VEP on the casing
- / analysis: **FFT, mapping, TPM, DSA/CSA, BSI** and others
- / meets the requirements of the **GDPR**
- / continuous development of the software
- / **backward compatibility**



ELMIKO BIOSIGNALS sp. z o. o.

Jeżewskiego 5c/7
02-796 Warsaw, Poland

+48 22 644 37 37
+48 600 37 37 90

elmiko@elmiko.pl
elmiko.pl



CREATE TO HELP

ISO13485 QUALITY MANAGEMENT SYSTEM AND CE MARK CERTIFIED BY TÜV RHEINLAND®